

DISPENSING INSTRUCTIONS

1. Select the frame. The frame should accommodate a minimum 13mm fitting height to the bottom of the eyewire and 10mm to the top. Adjust the frame for comfort and accuracy before taking measurements. Adjustable nose pads are recommended. Set the pantoscopic angle to 10-12°. Frame should have a slight face form.

2. PD and Fitting Height. Measure monocular fitting height by marking each demo lens at the pupil centers with a felt tip pen. Measure monocular PD using a pupillometer or by using the fitting height marks. To translate lens markings into measurements, place the frame on the center of the triangle, ensuring the marks on the lens are on the zero (0) line. Using the chart, record the monocular PD and monocular vertical heights.

3. Frame verification. Line up the pupillary mark on the demo lens with the cross on the chart. Verify that the distance and near zones are within the blue circle and that the eyewire is within the cutout diameter. This will ensure the minimum fitting height and cutout specifications are met. Confirm that the lens cutout is compatible with the material type you are specifying.

4. Include this information. Make sure you include the following information in your KODAK Unique Progressive lab order. These pieces of information are critical to the correct customization and optimization of the design to the frame size and shape.


- Monocular PD measurements
- Monocular fitting height measurements
- Manually traced **right** eyewire drawing
- Frame A, B, and DBL dimensions
- Frame brand, model, and eyesize

Note: KODAK Unique Progressives cannot be supplied unless items "a" through "d" above are provided.

5. Dispensing. Confirm the monocular PD and fitting height. Verify the lens R on the lensometer. Confirm the fit on the patient by verifying that the fitting cross is properly positioned over the pupil. Adjust the frame as necessary.


6. Teach proper viewing. Demonstrate the different viewing areas and appropriate head and eye movement.

PD and Fitting Height Chart




IMPORTANT If the entire lens is within the yellow area, it will be too small to be surfaced. You should recommend a different frame.

Frame Verification Chart (LEFT)



Frame Verification Chart (RIGHT)



Signet Armorlite, Inc. • 1001 Armorlite Dr. • San Marcos, CA 92069 • 800-950-5367 • www.signetarmorlite.com

Kodak
Licensed Product

Kodak and the Kodak trade dress are trademarks of Kodak, used under license by Signet Armorlite, Inc. SA PolarShades, CleAR and TLX 1.6 are trademarks and EvoClear, InstaShades and PolyClear are registered trademarks of Signet Armorlite, Inc. Trivex is a trademark of PPG Industries. Corning and SunSensors are registered trademarks of Corning Inc. Transitions is a registered trademark of Transitions Optical, Inc. Life Rx is a trademark of Vision-Ease Lens, Inc. Drivewear is a registered trademark of Younger Optics. ©2008 Signet Armorlite, Inc.